

## Varian GC Injector Liners

All items are subject to prior sale. Prices & availability on request !

old P/N	new P/N	Restek P/N	Description	Dimensions
<b>Liners for 1075 / 1077 Injectors</b>				
01-900109-00			empty	L = 72.0; ID = 2.0, OD = 6.3
01-900109-01		20792	with GlassWool	L = 72.0; ID = 4.0, OD = 6.3
01-900109-03		20715	with Frit	L = 72.0; ID = 4.0, OD = 6.3
01-900109-05		20721	Splitless	L = 74.0; ID = 2.0, OD = 6.3
01-900109-0X		20905-209.5	Splitless with CarbonFrit	L = 74.0; ID = 4.0, OD = 6.3
03-949809-00			with GlassWool	L = 72.0; ID = 4.0, OD = 6.3
03-949809-00			with GlassWool	L = 72.0; ID = 2.0, OD = 6.3
16-000830-00			Splitless, open insert	L = 72.0; ID = 4.0, OD = 6.3
16-000830-00			Splitless, open insert	L = 75.0; ID = 4.0; OD = 6.3
16-000830-01			with Frit	L = 72.0; ID = 4.0, OD = 6.3
03-949437-00			Splitless, open insert	L = 75.0; ID = 2.0; OD = 6.3
<b>Liners for 1078 / 1079 Injectors</b>				
01-900066-18			Splitless, empty	L = 54.0; ID = 2.00, OD = 5.0
03-918332-02			empty	L = 54.0; ID = 0.80, OD = 5.0
03-918332-03			with GlassWool	L = 54.0; ID = 0.80, OD = 5.0
03-918464-01	03-926119-46	21708	with Frit	L = 54.0; ID = 3.40, OD = 5.0
03-918466-00	03-926119-47	21711	Splitless, empty	L = 54.0; ID = 2.00, OD = 5.0
03-925330-00			empty	L = 54.0; ID = 0.XX, OD = 5.0
03-925331-00	03-926119-49	20992	empty	L = 54.0; ID = 0.50, OD = 5.0
03-925350-00			with GlassWool	L = 54.0; ID = 2.00, OD = 5.0
<b>Liners for SPI Injectors</b>				
01-900109-06		20775	Column ID 0.25 & 0.32	L = 54.0; ID = 0.53, OD = 4.6
01-900109-07		20778	Column ID 0.53	L = 54.0; ID = 0.80, OD = 4.6

01-900109-08

20850

Column ID 0.25, 0.32 & 0.53

L = 54.0; ID = 2.40, OD = 4.6